


25-09


Eye-tracking and pupillometry



When and why do we need eye-tracking?

Descriptive research: we're interested in the what, why and where of eye movements themselves

Explanatory research: oculomotor data may provide a window onto various cognitive processes



Terminology

Saccade



Saccadic amplitude

Saccadic latency

Fixation

Fixation duration

Microsaccade



Terminology

Saccade



Saccadic amplitude

Saccadic latency



Fixation

Fixation duration

Microsaccade




Eye position: two signals




Pupil location

Corneal reflection of (infrared) light sent from camera




Eye position: calibrate




Pupil location

Corneal reflection of (infrared) light sent from camera




Eye position: calibrate






Pupil location




Corneal reflection of (infrared) light sent from camera





In OpenSesame...





In OpenSesame...






In OpenSesame...




In OpenSesame...



In OpenSesame...




In OpenSesame...



Thus far:

Tracking the eye position (and pupil size) in a normal behavioral experiment

But what about **gaze-contingent** trickery?



From today's module on Canvas, download *eyetracking.oexp*

Mac users:
may have to install
PyGaze package
manually

Install from source

On other systems, you can install PyGaze as follows:


1. Download the PyGaze source code (.zip) from <https://github.com/esdalmaijer/PyGaze>.
 - Do *not* download the standalone Windows packages provided on the PyGaze website.
 - Verify that the version of PyGaze is compatible with your version of OpenSesame, as described [here](#).
2. Extract the .zip archive somewhere.
3. Inside, you will find these folders:
 - **opensesame_plugins**: As the name suggests, this folder contains the OpenSesame plugins, which need to be copied to (one of) the plugin folders, as described [here](#).
 - **pygaze**: This is the PyGaze Python library. You need to copy this to a folder in the Python path. On Windows, you can copy this folder to the OpenSesame program folder.
4. Done!

A variant of the Posner cueing task

Attention is biased by top-down cues


e.g., it takes longer for you to note the square on the right, when the arrow points to the left.

Task: indicate location of square (left / right)



Potential confound:

Simon effect (Maybe arrow biased response rather than attention)





A variant of the Posner cueing task

Attention is biased by top-down cues

e.g., it takes longer for you to note the square on the right, when the arrow points to the left.

Task: Move eyes to the square

Potential confound averted



A variant of the Posner cueing task

New experiment

experiment

instruction

pygaze_init

experiment_loop

trial_sequence

pygaze_drift_correct

pygaze_start_recording

stimulus


saccade_response

pygaze_stop_recording

feedback

logger

Unused items (0)



 **saccade_response** — inline script

Executes Python code

Prepare

Run

(080, 010)



```
1  # Here we determine the coordinates of the square, depending on condition:
2
3  - if exp.get('square') == 'left':
4      left_boundary = 192
5      right_boundary = 256
6  - else:
7      left_boundary = 768
8      right_boundary = 832
9
10 # And note that the upper and lower boundary are always the same (352 and 416).
```

A variant of the Posner cueing task

New experiment

experiment

instruction

pygaze_init

experiment_loop

trial_sequence

pygaze_drift_correct

pygaze_start_recording

stimulus

> saccade_response

pygaze_stop_recording

feedback

logger

Unused items (0)

saccade_response – inline script

Executes Python code

Prepare

Run

(033, 023)

```
1  # We want to know the response time, so first we have to mark the
2  # timepoint of stimulus onset:
3  start_time = self.time()
4
5  - while True: # Then we enter an endless loop, that we only break out of when
6    # the eyes land on the square
7
8    x,y = eyetracker.sample() # Check the eye position
9
10   # Check if the eye position is in the square:
11   - if x in range(left_boundary,right_boundary) and y in range(352,416):
12     # If so, we first want to make sure it's a true response, and not some
13     # accidental measurement (e.g. due to blinking). Let's say the response
14     # is real if the eyes are there for at least 50 ms. So we take a 50ms break.
15     self.sleep(50)
16     #... and then check again:
17     - if x in range(left_boundary,right_boundary) and y in range(352,416):
18       # If so, then specify the response time (i.e. current time minus start_time
19       # and minus the 50 ms pause)
20       exp.set('response_time',self.time() - start_time-50)
21       break # And break out of the loop
22   # Note that if the eye position isn't on the square anymore, we just start a
23   # new cycle of the loop.
```

Eye-tracker data

In our experiment we have a simple response_time variable...

Do we need more?

If possible, avoid
having to dive into
these files

→ define DV's
in OpenSesame


sub_7 - Notepad

File Edit Format View Help

```


ESACC L 2248653      2248710 58      627.1  534.1  1188.5  544.5  12.17  369
SFIx L  2248711
EFIX L  2248711      2249027 317     1186.3  548.4   1914
SSACC L 2249028
ESACC L 2249028      2249090 63      1185.3  548.4   703.4  522.3  10.48  357
SFIx L  2249091
MSG      2249659 move to word 2
EFIX L  2249091      2250585 1495    664.3  526.3   2158
SSACC L  2250586
SBLINK L 2250609
EBLINK L 2250609      2250686 78
ESACC L  2250586      2250742 157     667.2  524.4   698.0  543.8   0.79  1264
SFIx L  2250743
EFIX L  2250743      2251010 268     694.0  544.1   2347
SSACC L  2251011
ESACC L  2251011      2251055 45      692.8  548.2  1110.7  568.5   9.10  421
SFIx L  2251056
EFIX L  2251056      2251299 244     1124.1  573.1   1828
SSACC L  2251300
ESACC L  2251300      2251340 41      1129.8  571.0  1421.5  569.2   6.19  304
SFIx L  2251341
EFIX L  2251341      2251749 409     1429.4  569.5   1515
SSACC L  2251750
ESACC L  2251750      2251821 72      1432.3  565.8   700.4  535.4  15.64  418
SFIx L  2251822
MSG      2252601 move to word 3
EFIX L  2251822      2253059 1238    675.0  533.8   1626
SSACC L  2253060
ESACC L  2253060      2253105 46      670.2  542.0  1097.7  566.2   9.32  378
SFIx L  2253106
EFIX L  2253106      2253207 202     1105.1  560.8   1000

```



Eye-tracker data

→ define DV's *in* OpenSesame
Example: saccadic curvature




General properties | saccade_response

saccade_response — inline script
Executes Python code


Prepare Run (001, 001)

```
1 # We want to know the response time, so first we have to mark the
2 # timepoint of stimulus onset:
3 start_time = self.time()
4
5 all_y_values = []
6
7 while True: # Then we enter an endless loop, that we only break out of when
8             # the eyes land on the square
9
10     x,y = eyetracker.sample() # Check the eye position
11
12     all_y_values.append(384-y) # add deviation from vertical center to list
13
14     # Check if the eye position is in the square:
15     if x in range(left_boundary,right_boundary) and y in range(352,416):
16         # If so, we first want to make sure it's a true response, and not some
17         # accidental measurement (e.g. due to blinking). Let's say the response
18         # is real if the eyes are there for at least 50 ms. So we take a 50ms break.
19         self.sleep(50)
20         #... and then check again:
21         if x in range(left_boundary,right_boundary) and y in range(352,416):
```



Eye-tracker data

→ define DV's *in* OpenSesame
Example: saccadic curvature




```
General properties | saccade_response
```

saccade_response – inline script
Executes Python code

Prepare Run (001, 001)


```
1 # We want to know the response time, so first we have to mark the
2 # timepoint of stimulus onset:
3 start_time = self.time()
4
5 all_y_values = []
6
7 - while True: # Then we enter an endless loop, that we only break out of when
8   # the eyes land on the square
9
10    x,y = eyetracker.sample() # Check the eye position
11
12    all_y_values.append(384-y) # add deviation from vertical center to list
13
14    # Check if the eye position is in the square:
15    - if x in range(left_boundary,right_boundary) and y in range(352,416):
16      # If so, we first want to make sure it's a true response, and not some
17      # accidental measurement (e.g. due to blinking). Let's say the response
18      # is real if the eyes are there for at least 50 ms. So we take a 50ms break.
19      self.sleep(50)
20      #... and then check again:
21      - if x in range(left_boundary,right_boundary) and y in range(352,416):
```




Eye-tracker data

→ define DV's *in* OpenSesame
Example: saccadic curvature

```
12 all_y_values.append(384-y) # add deviation from vertical center to list
13
14 # Check if the eye position is in the square:
15 - if x in range(left_boundary,right_boundary) and y in range(352,416):
16     # If so, we first want to make sure it's a true response, and not some
17     # accidental measurement (e.g. due to blinking). Let's say the response
18     # is real if the eyes are there for at least 50 ms. So we take a 50ms break.
19     self.sleep(50)
20     #... and then check again:
21 - if x in range(left_boundary,right_boundary) and y in range(352,416):
22     # If so, then specify the response time (i.e. current time minus start_time
23     # and minus the 50 ms pause)
24     exp.set('response_time',self.time() - start_time-50)
25     break # And break out of the loop
26     # Note that if the eye position isn't on the square anymore, we just start a
27     # new cycle of the loop.
28
29
30 curvature = int(float(sum(all_y_values))/len(all_y_values))
31 exp.set('curvature',curvature)
32
```



Pupillometry practical



Pupillometry practical

Assignment due: Sunday 1st 23:59

Pupillary light response: not just to our direct visual environment



Also to memorized brightness

(Mathôt et al., 2017: pupil response to semantic brightness of words)



Pupillometry practical


Assignment due: Sunday 1st 23:59




Pupillometry practical

Assignment due: Sunday 1st 23:59

“have you seen this
orientation?”



H: pupil responds to brightness
of memorized stimulus location



Pupillometry practical

Assignment due: Sunday 1st 23:59

*From today's module in Canvas, download
pupillometry practical.pdf*

Assignment consists of 2 parts; Part 2
is about analyzing data (*Thursday!*)

